GeneChip® Human Transcriptome Array 2.0

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Controls	
ERCC probe set ^{1,2}	63
Background probes	Antigenomic set
Poly-A controls ²	<i>dap, lys, phe, thr</i>
Hybridization controls	<i>bioB, bioc, bioD, creX</i>

¹ Probe sets interrogating external RNA controls present in the Ambion® ERCC RNA Spike-In Control Mixes, P/N 4456740 and 445673.

² This array contains probe sets for both ERCC and Poly-A spike in-controls. Sequence homology between the two control mixes will result in cross-hybridization of target to the control probes on the array. It is important to use only one control probe set when processing the array (ERCC or Poly-A controls) but not both.

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货号	产品	描述
902233	GeneChip® Human Transcriptome Array 2.0	Contains 2 arrays
902162	GeneChip® Human Transcriptome Array 2.0	Contains 10 arrays

Related products

货号	产品	描述
901524	GeneChip® WT Terminal Labeling and Controls Kit ¹	Sufficient for 30 reactions
901525		Sufficient for 10 reactions
900670	GeneChip® WT Terminal Labeling Kit	Sufficient for 10 reactions
900671		Sufficient for 30 reactions
900433	GeneChip® Poly-A Control Kit	Sufficient for 100 reactions
900454	GeneChip® Hybridization Control Kit	Sufficient for 30 reactions
900720	GeneChip® Hybridization, Wash, and Stain Kit	Sufficient for 30 reactions

¹ Contains Terminal Labeling Kit, Poly-A Control Kit, and Hybridization Control Kit

www.affymetrix.com Please visit our website for international distributor contact information.

“For Research Use Only. Not for use in diagnostic procedures.”

P/N EXP01883 Rev. 3

©Affymetrix, Inc. All rights reserved. Affymetrix®, Axiom®, Command Console®, CytoScan®, DMET™, GeneAtlas®, GeneChip®, GeneChip-compatible™, GeneTitan®, Genotyping Console™, myDesign™, NetAffx®, OncoScan™, Powered by Affymetrix™, PrimeView®, Procarta®, and QuantiGene® are trademarks or registered trademarks of Affymetrix, Inc. All other trademarks are the property of their respective owners.

Products may be covered by one or more of the following patents: U.S. Patent Nos. 5,445,934; 5,744,305; 5,945,334; 6,140,044; 6,399,365; 6,420,169; 6,551,817; 6,733,977; 7,629,164; 7,790,389 and D430,024 and other U.S. or foreign patents. Products are manufactured and sold under license from OGT under 5,700,637 and 6,054,270.

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